


Exome Sequencing in genetic diseases


The clinical / clinical pathology partnership that benefits patients

Exome Sequencing in diagnostic practice

The term “genetic disease” embraces a wide range of pathologies. The identification of causal variants facilitates clinical diagnosis and prognosis, patient management, genetic counselling and, in some cases, the adjustment of therapeutic care to individual patient needs.



The conventional step-by-step sequencing approach (gene by gene or panel by panel), which sometimes follows complex decision trees, is gradually being replaced by **exome sequencing as a first-tier technique, in an increasing number of indications**¹⁻².



Exome Sequencing makes it possible to sequence, in a single step, all of the coding regions of the human genome (containing more than 85% of pathogenic mutations) **at an affordable price.**

As a first-tier technique, exome sequencing provides:



time-saving compared to carrying out a series of gene panels,



a diagnostic test for patients where the clinical examination does not enable the clinician to identify the gene/panel of genes to be tested.



retrophenotyping in atypical clinical pictures,



an improved diagnostic yield compared to the gene panel.

The exome sequencing by Eurofins Biomnis

4 weeks

In the laboratory

- A dedicated team of clinical pathologists, technicians and bioinformatics specialists
- Continuous biopathological support, from the implementation of the test to the interpretation of the results.

Sequencing

- Whole exome coverage: ~22,000 genes and 37.5 Megabases targeted
- Depth > 30X for ~98% bases*
- >99% of recall**

At the end of the analytical process

- An interpretation of an average of 30,000 variants in consultation with the clinician
- Provision of raw data (fastQ, VCF, BAM and quality report) via a secure interface
- Report with detailed results within an optimised period of 4 weeks

Why Eurofins Biomnis?



- Specialised medical pathology laboratory
- ISO 15189 accreditation (CAP standard equivalent)
- Authorisation to perform constitutional genetics tests
- Certified clinical pathologists
- Expertise in sequencing techniques in diagnostic practice

The clinical/clinical pathology partnership that benefits patients

One of the challenges of exome sequencing is the interpretation of variants.


To actively involve partner clinicians and clinical pathologists in the interpretation of data, Eurofins Biomnis provides access to **SeqOne**, a secure bioinformatics platform.

*target CDS Refseq +/- 3 base pairs

**data calculated from SNV's from NIST002 reference samples, for 40 million pairs of reads generated.

SeqOne, for optimal diagnostic performance

Thanks to this platform, a **joint interpretation with Eurofins Biomnis or an autonomous interpretation** is possible. This ensures **optimal diagnostic performance** (laboratory experience and literature³).



ACMG	VEB	Gene	Transcript	NA Mutation	AA Mutation	Impact	Effect	Transmission	Databases	Pop. freq.	Sample Freq. Project
Pathogenic	Pathogenic	SLC25A22	NM_024096.5	c.3712_3714del	p.A1770V (Ter Stop)	High	stopgain variant	[AF]	dbSNP	1%	2/3 66.67%
Pathogenic	Pathogenic	SLC25A22	NM_024096.5	c.3710G>A	p.R127K	Missense	missense variant	[AF]	dbSNP	41 x 10 ⁻⁴	2/3 66.67%
Pathogenic	Pathogenic	GGX	NM_002388.3	c.1057T>G		High	splice donor variant	[AF]	dbSNP	1.1 x 10 ⁻³	2/3 66.67%
Pathogenic	Pathogenic	WDR91	NM_024096.3	c.1704A>G	p.R171R	Missense	missense variant	[AF]	dbSNP	1%	2/3 66.67%
Pathogenic	Pathogenic	PIB	NM_002093.2	c.147C>G	p.P16R	Missense	missense variant	[AF]	dbSNP	1%	1/3 33.33%
Pathogenic	Pathogenic	STT3	NM_010284.3	c.482G>A	p.K128I	Missense	missense variant	[AF]	dbSNP	1.6 x 10 ⁻⁴	2/3 66.67%
Pathogenic	Pathogenic	ALM5	NM_010235.4	c.1023T>C	p.D105P	Missense	missense variant	[AF]	dbSNP	1%	2/3 66.67%
Pathogenic	Pathogenic	PCN1	NM_004281.5	c.2037C>G	p.R101G	Missense	missense variant	[AF]	dbSNP	1.2 x 10 ⁻³	2/3 66.67%
Pathogenic	Pathogenic	WDR91	NM_024096.3	c.1611_1612del		High	splice donor variant	[AF]	dbSNP	1%	1/3 33.33%
Pathogenic	Pathogenic	CFP8	NM_010284.3	c.228C>T	p.R167W	High	stop gained	1%	dbSNP	1.6 x 10 ⁻⁴	2/3 66.67%
Pathogenic	Pathogenic	ABMF12	NM_010284.4	c.368G>T	p.T112I	Missense	missense variant	[AF]	dbSNP	1.6 x 10 ⁻⁴	1/3 33.33%

The data is accessible and available to the clinician at all times, which allows:

- reanalysis, when if required;
- use of data for research.


Eurofins Biomnis exome add-on modules:

- **CNV detection**
- **Detection of mitochondrial variants**, even with low heteroplasmy
- **Detection of pathogenic variants** in the *MUC1* VNTR
- **Detection of uniparental disomy** (for trio exomes)
- **Pharmacogenetic analysis**
- **HLA typing**



Sequence your exomes with Eurofins Biomnis

3 interpretation options/solutions



References

1. Exome sequencing has higher diagnostic yield compared to simulated disease-specific panels in children with suspected monogenic disorders. Oliver James Dillon et al. European Society of Human Genetics, Feb 18, 2018.
2. Meta-analysis and multidisciplinary consensus statement: exome sequencing is a first-tier clinical diagnostic test for

individuals with neurodevelopmental disorders. Srivastava S et al. Genet Med. 2019 Nov;21(11):2413-2421. doi: 10.1038/s41436-019-0554-6. Epub 2019 Jun 11.

3. Paediatric genomics: diagnosing rare disease in children. Wright CF et al. Nat Rev Genet. 2018 May;19(5):253-268. doi: 10.1038/nrg.2017.116. Epub 2018 Feb 5.



Practical details

	Exome Sequencing & Access to the interpretation interface	
Interpretation level	Detailed report	Basic report
Turnaround	<ul style="list-style-type: none"> ● 4 weeks excluding any additional examinations 	2 weeks
Indications	<ul style="list-style-type: none"> ● intellectual disability, neurodevelopmental disorders, ● syndrome-based disorders, ● clinical pictures which do not directly suggest a specific gene test or gene panel, or negative result for these primary tests, ● organ damage (kidney, heart, etc.). 	
Sample	Solo: 5 mL EDTA whole blood or DNA sample Duo/Trio: 5mL EDTA whole blood or DNA sample from relative	
Conservation & transport	Room temperature	
Technique	Exome sequencing + bioinformatics pipeline	
Required documents	B34-INTGB Test request form available on www.eurofins-biomnis.com > Test guide > Analysis code EXOME	
Price	Contact us	
Complementary test	Study of relative persons by Sanger sequencing: 5 mL EDTA whole blood or DNA sample from relative	

For more information:

Eurofins Biomnis International Division
 17/19 avenue Tony Garnier
 BP 7322 - 69357 LYON Cedex 07 - FRANCE
 E-mail: international@biomnis.eurofinseu.com
www.eurofins-biomnis.com



Biomnis